

NK Cell Degranulation Assay

Reagents

Reagent	Vendor	Catalogue #	Stock Conc.
FACS Tubes	BD Falcon	352054	-
P815 Ab=rabbit anti-mouse lymphocyte	Accurate Chemical	AIA3940	-
GolgiPlug (Brefeldin A (BFA))	BD	555029	-
GolgiStop (Monensin)	BD	554724	-
Cytofix/Cytoperm Kit	BD	554722	-
PMA	Sigma	1585	1mg/mL
Ionomycin	Sigma	141128	1mg/mL
RPMI-1640	Sigma	R0883	-
Pen/Strep	Cellgro	30-001-C1	Pen: 5000 IU/mL Strep: 5000ug/mL
L-glutamine	Cellgro	25-002-C1	200mM; 29.2mg/mL
HEPES (1M; 238.3mg/mL)	Cellgro	25-060-C1	1M
FBS, Heat-inactivated	Sigma	F4135	-
PBS	Sigma	D8537	-

Note: addition of GolgiStop is necessary AND sufficient to look at CD107a as a marker for degranulation. A combination of GolgiStop and BFA improves intracellular staining for cytokines such as IFN-gamma. Please keep in mind that CD69 expression is blocked by BFA, and that expression of some surface markers is affected by these inhibitors.

Media:

R10

Reagent	Stock Concentration	Volume to Add	Final Concentration
RPMI 1640	-	500	-
FBS	100% (v/v)	55ml	10% (v/v)
Penicillin/Strep	Pen: 5000 IU/mL Strep: 5000ug/mL	5.5mL	Pen: 50 IU/mL Strep: 50ug/mL
L-glutamine	200mM	5.5mL	2mM
HEPES	1M	5.5ml	10mM

PMA

Working solution: 1mg/mL. Reconstitute in DMSO (1mg in 1mL)
 Prepare 20uL aliquots and keep at -20C
 Final concentration: 100ng/mL. Prepare a 1:10 dilution and add 1uL/mL culture

Ionomycin

Working solution: 1mg/mL. Reconstitute in DMSO (1mg in 1mL)
 Prepare 10uL aliquots and keep at -20C

Final concentration: 1ug/mL. Add 1uL/mL culture

Assay

Note: *The assay should ideally be performed using fresh PBMCs that have been isolated within 6h following collection in ACD tubes. However frozen PBMCs can be used.*

1. Isolate PBMCs as described in SOP #08-00

Note: *If necessary, this assay can be performed using enriched PBMC subpopulations such as CD3-depleted PBMCs or NK cells. In this case, use the relevant RosetteSep cocktail from StemCell Technologies according to manufacturer instructions before the ficoll procedure.*

2. Count PBMCs and resuspend the pellet in R10 at 1×10^6 /mL.

3. Count target cells from stock in culture

Note: *Here is a list of target cells commonly used to test NK cell function. These can be substituted by any other relevant cell line.*

K562 cells: *MHC-I-deficient, express NKG2D ligands, maintained in R10*

221 cells: *MHC-I-deficient, express NCR ligands, maintained in R10*

P815 cells: *murine cell line expressing Fc receptors, can be either coated with p815 antibodies in order to assess ADCC function, or coated with any specific Ab to perform a redirected-lysis assay. To coat p815 cells, add 10ug/mL Ab (10uL of stock) to 1 million p815 cells in 1 mL R10, incubate for 1h at 37C, wash and resuspend pellet in 1mL R10. In parallel, prepare 1 million of uncoated p815 to use as control.*

4. Put 1 million target cells in a FACS tube or a 15mL conical (adjust amount of cells if more than 10 patients)

5. Spin 350xg, 10min. Discard supernatant

6. Resuspend pellet at 1×10^6 /1mL R10

7. Put 1mL PBMCs in each reaction tube (FACS tube) and set up your experiment as follows:

- a. Unstimulated control: no target cells
- b. Stimulation control: add 1uL of PMA and 1uL of ionomycin
- c. Add 100uL of 221 cells
- d. Add 100uL of K562 cells
- e. Add 100uL of uncoated p815
- f. Add 100uL of Ab-coated p815 cells
- g. Any other cell line to be tested.....

Note: *100uL of target cells results in an E:T ratio of 10:1 which is optimal when using whole PBMCs and the abovementioned target cells. If testing other stimuli or cell lines, it is recommended to try other ratios (5:1 or 20:1). If using pure NK cells, start with 1:1 or 2:1 E:T ratios.*

8. In each tube, add 1uL of GolgiStop/mL and a CD107a conjugated antibody. If looking at cytokine-producing cells by ICS, add 1uL BFA/mL. If more than 3 tubes, prepare a master mix.

9. Incubate the PMA-ionomycin reaction tube for 2h at 37C, then put at 4C protected from light. Longer incubation time results in massive cell death.

10. Incubate other reaction tubes for 4-6H at 37C, then proceed with the staining or put at 4C O/N protected from light. **If stimulation requires more than 16h of incubation, add the CD107a Ab and GolgiStop/BFA only for the last 6h of incubation.**